Instructions for Use

IMPORTANT INFORMATION — PLEASE READ

Caution: U.S. federal law restricts this device to sale by, or on the order of, a licensed dentist or physician.

■ General Information

The Inclusive Tapered Implant System consists of dental implants, prosthetic components, surgical instrumentation, and related accessories packaged under the Inclusive brand name for use by qualified, licensed clinicians and laboratory technicians fully trained in their application.

For specific product identification and contents, please refer to individual product labels and the following catalog:

- Inclusive Dental Implant System Product Catalog (**MKT 787**)

For detailed information on the specifications and intended use of a particular product, please refer to the following user manuals:

- Inclusive Dental Implant System Surgical Manual (**UM 4236**)
- Inclusive Prosthetic Components Restorative Manual (**UM 5088**)

■ Explanation of Label Symbols

The symbols glossary is provided on page 9 of this IFU document.

■ Disclaimer of Liability

The guidelines presented herein are not adequate to allow inexperienced clinicians to administer professional implant treatment or prosthetic dentistry, and are not intended to substitute for formal clinical or laboratory training. Inclusive devices should only be used by individuals with training and experience specific to their clinically accepted application.

Prismatik Dentalcraft, Inc. is not liable for damages resulting from treatment outside of our control. The responsibility rests with the provider.

■ MRI

The Inclusive Tapered Implant System has not been evaluated for safety and compatibility in the magnetic resonance (MR) environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the Inclusive Tapered Implant System in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

DENTAL IMPLANTS

■ Description

Inclusive Tapered Implants are endosseous devices manufactured from titanium alloy. They are compatible with the prosthetic components and surgical instrumentation of the Inclusive Tapered Implant System.

■ Indications for Use

Inclusive Tapered Implants are indicated for use in maxillary and mandibular partially or fully edentulous cases, to support single, multiple-unit, and overdenture restorations. The implants are to be used for immediate loading only in the presence of primary stability and appropriate occlusal loading.

■ Contraindications

Inclusive Tapered Implants should not be placed in patients discovered to be medically unfit for the intended treatment. Prior to clinical intervention, prospective patients must be thoroughly evaluated for all known risk factors and conditions related to oral surgical procedures and subsequent healing. Contraindications include but are not limited to:

- vascular conditions
- uncontrolled diabetes
- clotting disorders
- anticoagulant therapy
- metabolic bone disease
- chemotherapy or radiation therapy
- chronic periodontal inflammation
- insufficient soft tissue coverage
- metabolic or systemic disorders associated with wound and/or bone healing
- use of pharmaceuticals that inhibit or alter natural bone remodeling
- any disorders which inhibit a patient's ability to maintain adequate daily oral hygiene
- uncontrolled parafunctional habits
- insufficient height and/or width of bone, and insufficient interarch space.

Treatment of children is not recommended until growth is finished and epiphyseal closure has occurred.

■ Warnings

- Do not reuse Inclusive Tapered Implants. The reuse of such device on another patient is not recommended due to the risks of cross-contamination or infection.
- Inclusive Tapered Implants may only be used for their intended purpose in accordance with general rules for dental/surgical treatment, occupational safety, and accident prevention. They must only be used for dental procedures with the restorative components they were designed for. If the indications and intended use are not clearly specified, treatment should be suspended until these considerations have been clarified.

- The following instructions are not sufficient to allow inexperienced clinicians to administer professional prosthetic dentistry. Inclusive Tapered Implants, surgical instruments, and restorative components must only be used by dentists and surgeons with training/experience with oral surgery, prosthetics and biomechanical requirements, as well as diagnosis and preoperative planning.
- The implant site should be inspected for adequate bone by radiographs, palpations and visual examination. Determine the location of nerves and other vital structures and their proximity to the implant site before any drilling to avoid potential injury, such as permanent numbness to the lower lip and chin.
- Absolute success cannot be guaranteed. Factors such as infection, disease, and inadequate bone quality and/or quantity can result in osseointegration failures following surgery or initial osseointegration.

■ Precautions

Surgical Procedures

Minimizing tissue damage is crucial to successful implant osseointegration. In particular, care should be taken to eliminate sources of infection, contaminants, surgical and thermal trauma. Risk of osseointegration failure increases as tissue trauma increases. For best results, please observe the following precautions:

- All drilling procedures should be performed at 2000 RPM or less under continual, copious irrigation.
- All surgical instruments used must be in good condition and should be used carefully to avoid damage to implants or other components. Implants should be placed with sufficient stability; however, excessive insertion torque may result in implant fracture, or fracture or necrosis of the implant site. The proper surgical protocol should be strictly adhered to.
- Since implant components and their instruments are very small, precautions should be taken to ensure that they are not swallowed or aspirated by the patient.
- Prior to surgery, ensure that the needed components, instruments and ancillary materials are complete, functional and available in the correct quantities.

Prosthetic Procedures

Following successful placement of Inclusive Tapered Implants, verify primary stability and appropriate occlusal loading before proceeding with the placement of a permanent or provisional prosthesis. All components that are used intraorally should be secured to prevent aspiration or swallowing. Distribution of stress is an important consideration. Care should be taken to avoid excessive loads significantly transverse to the implant axes.

■ Sterility

Inclusive Tapered Implants are shipped sterile. They should not be resterilized. They are for single use only, prior to the expiration date. Do not use implants if the packaging has been compromised or previously opened.

■ Storage and Handling

Inclusive Tapered Implants must be stored in a dry location at room temperature, in their original packaging. Inclusive Tapered Implants are packaged sterile. Do not handle implant surfaces directly. Users are advised to visually inspect packaging to ensure seals and contents are intact prior to use. Please refer to the individual product label for all relevant product information and cautions.

■ INSTRUCTIONS FOR USE

Soft Tissue Reflection

Following administration of anesthesia, make an incision designed for elevation of a flap. Perform alveoloplasty on the crest of the ridge, if needed, to create a more even plane in which to place the implant. Irrigation should be used for all modifications of the bone.

Site Preparation

Step 1: Lance Drill – With copious irrigation, perforate the alveolar crest. Utilize a surgical guide, if necessary, as a reference for proper positioning.

Step 2: Pilot Drill – Select the 2.3/2.0 mm Pilot Drill. If any change is needed in trajectory, it may be corrected at this time. With copious irrigation, drill a pilot hole to the appropriate depth marking on the drill.

Check the orientation of the initial osteotomy using a Parallel Pin. If placing more than one implant and parallelism is desired, begin drilling the next site and align as the trajectory of the bone permits.

Step 3: Surgical Drills – Depending on implant diameter and the density of bone at the osteotomy site, it may be necessary to utilize one or more of the Surgical Drills to widen the osteotomy. To avoid over-preparation, widening drill diameters should be used only as necessary, and in proper succession. Select the desired Surgical Drill, accounting for the density of bone at the osteotomy site and the diameter of the implant to be placed. With copious irrigation, drill to the appropriate depth marking on the drill. The final drill for each implant diameter should be based on bone density (soft or hard), as charted below. The goal is to achieve high primary stability upon implant placement.

Drilling Sequence Chart							
Implant Size	Lance Drill (Ø1.5 mm)	Pilot Drill (Ø2.3/2.0 mm)	Surgical Drill (Ø2.8/2.3 mm)	Surgical Drill (Ø3.4/2.8 mm)	Surgical Drill (Ø3.8/3.4 mm)	Surgical Drill (Ø4.4/3.8 mm)	Surgical Drill (Ø4.9/4.4 mm)
Ø3.2 mm	Soft/Hard	Soft/Hard	Hard Only				
Ø3.7 mm	Soft/Hard	Soft/Hard	Soft/Hard	Hard Only			
Ø4.2 mm	Soft/Hard	Soft/Hard	Soft/Hard	Soft/Hard	Hard Only		
Ø4.7 mm	Soft/Hard	Soft/Hard	Soft/Hard	Soft/Hard	Soft/Hard	Hard Only	
Ø5.2 mm	Soft/Hard	Soft/Hard	Soft/Hard	Soft/Hard	Soft/Hard	Soft/Hard	Hard Only

Step 4: (Optional) Dense Bone Tap – If indicated by the presence of dense bone, select the Screw Tap with a diameter matching that of the implant. Place the tap into the prepared implant site. Apply firm pressure and begin slowly rotating the tap (25 RPM maximum). When the threads begin engaging the bone, allow the tap to feed into the site without applying additional pressure. The osteotomy should be tapped through the cortical bone. Reverse the tap out of the site.

Implant Placement

Step 1: Implant Selection – Remove the titanium implant holder from its packaging and place it onto a sterile field.

Step 2: Initial Placement – Use slight finger pressure to pinch the occlusal end of the implant in its holder while inserting the appropriate Implant Driver. Gently rotate implant and holder, allowing the driver to engage the implant connection. With the driver securely attached to the implant, squeeze the opposing end of the holder to disengage the implant from the holder. Transport the implant to the prepared site, and insert into the osteotomy. Rotate clockwise with applied pressure to engage the self-tapping grooves. Avoid lateral forces, which can affect the angulation and final alignment of the implant.

Step 3: Advancement and Final Seating – Continue threading the implant into the osteotomy site using the preferred placement method. A minimum torque value of 35 Ncm upon final seating indicates good primary stability.

Methods of Implant Placement

Option 1: Handpiece Implant Placement – Place the appropriate Implant Driver into the handpiece. Seat the driver into the internal hex connection of the implant, and press firmly to fully engage the connection. Thread the implant into the osteotomy at approximately 25 RPM until fully seated.

Option 2: Manual Implant Placement – Assemble the Adjustable Torque Wrench with the Surgical Adaptor and appropriate Implant Driver. With the implant threaded securely in its site, seat the driver into the internal hex connection of the implant, and press firmly to fully engage the connection. Turn the wrench clockwise in increments of approximately 90 degrees. Avoid lateral forces, which can affect final alignment of the implant.

Implant Positioning

The implant should be rotated at the time of placement to ensure optimal positioning of the internal hex connection. This will allow the restoring clinician to take full advantage of the anatomical abutment contours and minimize the need for abutment preparation. Adjust the final position of the implant so that any one of the six flats of the internal hex connection is oriented toward the facial.

Bone Profiling

If indicated by subcrestal placement of the implant or excess bone around the restorative platform, select the appropriate Bone Profiler based on the implant's platform size. Insert the tip of the profiler into the implant's connection interface, taking care to ensure a parallel orientation with the implant. Using finger pressure, maintain parallelism while gently rotating the profiler clockwise until the profiler reaches a stop against the implant platform (progressing no farther with continued rotation). Remove the profiler when finished. With suction, irrigate site to ensure removal of bone debris.

NOTE: Failing to maintain parallelism between the profiler and the implant may result in damage to the implant interface, or lead to incorrect profiling of the bone.

Healing Component Placement

Following implant placement, prepare the site for healing by placing either a Healing Abutment (single-stage surgical protocol) or the Cover Screw (two-stage surgical protocol).

Option 1: Healing Abutment – If observing a single-stage surgical protocol, select a Healing Abutment of the appropriate height and diameter. Thread the abutment into place atop the implant. Hand-tighten with the appropriate Prosthetic Driver.

Option 2: Cover Screw – If observing a two-stage surgical protocol, thread the Cover Screw into place atop the implant. Hand-tighten with the appropriate Prosthetic Driver.

Closure and Suturing

If the soft tissue was reflected, close and suture the flap utilizing the desired technique. Take a postoperative radiograph to use as a baseline, and advise the patient as to the recommended postoperative procedures.

Second-Stage Uncovery (Two-Stage Surgical Protocol)

Following the appropriate healing period, make a small incision in the gingiva over the implant site to expose the Cover Screw. Using the Prosthetic Driver, remove the Cover Screw, and place a Healing Abutment or Temporary Abutment of the appropriate height and diameter.

PROSTHETIC COMPONENTS

■ Description

Prosthetic components of the Inclusive Tapered Implant System, consisting of abutments, screws, analogs, copings, and related restorative accessories, are manufactured from titanium alloy, gold alloy, polymers, or zirconia. They are compatible with Inclusive Tapered Implants. For product-specific descriptions and sterility information, please refer to the individual product labels and appropriate Inclusive catalog and/or user manuals.

■ Indications for Use

Inclusive Abutments are premanufactured prosthetic components directly connected to endosseous dental implants and are intended for use as an aid in prosthetic rehabilitation.

All digitally designed abutments for use with Inclusive Abutments for CAD/CAM are intended to be manufactured at a validated milling center.

Inclusive Multi-Unit Abutments are intended to provide support and retention for multi-unit screw-retained restorations. The 30-degree multi-unit abutments must be used within 45 degrees of parallelism for a splinted restoration. The 17-degree multi-unit abutments must be used within 32 degrees of parallelism for a splinted restoration.

■ Contraindications

Inclusive Abutments

- Wall thickness less than 0.5 mm
- Gingival margin diameter less than 0.5 mm wider than implant
- Gingival height greater than 6 mm**
- Angle corrections of more than 30 degrees
- Angle corrections of more than 20 degrees*
- Less than 0.5 mm margin height
- Less than 4.0 mm abutment post height above the gingival collar

*ONLY applicable to Dentsply Implants Astra Tech Implant System® EV, Dentsply Implants Astra Tech OsseoSpeed™ and Camlog® Screw-Line.

**ONLY applicable to Dentsply Implants Astra Tech Implant System® EV.

Inclusive Hybrid Abutments (titanium bases with zirconia copings)

- Wall thickness less than 0.5 mm
- Gingival margin diameter less than 0.5 mm wider than implant
- Gingival height greater than 6.0 mm

- Angle corrections
- Less than 0.5 mm margin height
- Less than 4.0 mm abutment post height above the gingival collar

NOTE: A patient-specific finished device will consist of both the titanium base and zirconia coping. The titanium base is not intended to provide angle correction. The zirconia coping is intended to be milled straight only.

Inclusive Titanium Esthetic Abutments

- Wall thickness less than 0.5 mm
- Gingival margin diameter less than 0.5 mm wider than implant
- Angle corrections of more than 15 degrees
- Less than 0.5 mm margin height
- Less than 4.0 mm abutment post height above the gingival collar

NOTE: Small diameter implants with angled abutments are not recommended for the posterior region of the mouth.

Inclusive Multi-Unit Abutments

- Greater than 45 degrees divergence from parallel for a splinted restoration when using 30-degree multi-unit abutments
- Greater than 32 degrees divergence from parallel for a splinted restoration when using 17-degree multi-unit abutments

■ **Warnings**

An Inclusive abutment is intended to be used on an individual patient only. The reuse of such device on another patient is not recommended due to the risks of cross-contamination or infection. Small-diameter implants with angled abutments are not recommended for the posterior region of the mouth.

■ **Adverse Effects**

The following adverse effects have been observed when using prosthetic components and accessories:

- Components used in the patient's mouth have been aspirated or swallowed.
- The abutment screw has fractured due to application of excessive torque.
- The abutment is not adequately secured due to inadequate application of torque.

■ **Precautions**

Inclusive abutments may only be used for their intended purpose in accordance with general rules for dental/prosthetic treatment, occupational safety, and accident prevention. Inclusive abutments must only be used for dental procedures with the implant systems they were designed for. If the indications and intended use are not clearly specified, treatment should be suspended until these considerations have been clarified. All components that are used intraorally must be secured to prevent aspiration or swallowing. Prior to placement, ensure that the required components, instruments, and ancillary materials are complete, functional, and available in the correct quantities.

■ **Side Effects**

No side effects, according to current knowledge.

■ **Storage and Handling**

Inclusive Prosthetic Components labeled STERILE should be stored in a dry location at room temperature, in their original packaging. Users are advised to visually inspect packaging to ensure seals and contents are intact prior to use. Please refer to the individual product label for all relevant product information and cautions. Sterile products are intended for single-use only, prior to the expiration date. Do not use sterile products if the packaging has been compromised or previously opened. Do not resterilize.

Products labeled NON-STERILE should be cleaned and sterilized according to a validated method prior to use in the oral environment.

■ **Sterility**

Healing abutments and multi-unit abutments are shipped sterile. They should not be resterilized. They are for single use only, prior to the expiration date.

Non-sterile abutments and screws must be cleaned and sterilized prior to clinical use, according to a validated method.

- **Cleaning:** Prepare cleaning solution using 5 mL of dish soap per gallon of tap water. Fully immerse the devices in solution and scrub them with a soft-bristle brush. Remove the components and rinse them under running tap water. Dry the devices with a clean, lint-free cloth.

The recommended sterilization process is based on the ANSI/AAMI/ISO 17665-1 and ANSI/AAMI ST79 guidelines, as follows:

- **Sterilization of Inclusive Titanium and Inclusive Esthetic Abutments:** Gravity-fed sterilizers: Autoclave in sterilization pouch for 30 minutes at 121°C (250°F). Devices are to be used immediately after sterilization.
- **Sterilization of Inclusive Hybrid Abutments:** Gravity-fed sterilizers: Autoclave in sterilization pouch for 15 minutes at 132°C (270°F). Allow sterilized components to dry for at least 30 minutes.

NOTE: The validated procedures require the use of FDA-cleared sterilization trays, wraps, biological indicators, chemical indicators, and other sterilization accessories labeled for the sterilization cycle recommended. The healthcare facility should monitor the sterilizer for the facility according to an FDA-recognized sterility assurance standard such as ANSI/AAMI ST79.

■ **Dental Implant Compatibility**

Prosthetic components for the Inclusive Tapered Implant System are compatible with Inclusive Tapered Implants. The platform-specific compatibility of each component is indicated on the individual product label. The availability of a particular type of prosthetic component may be limited by restorative platform, geographical territory, or other considerations. For a complete product listing, please refer to the *Inclusive Dental Implant System Product Catalog*, or contact an Inclusive sales representative.

■ Recommended Torque Values

The recommended torque value for affixing Inclusive Tapered Implant abutments and multi-unit abutments to Inclusive Tapered Implants is indicated in the following table:

Platform Size	Manufacturer's Recommended Torque (Ncm)		
	Titanium Abutment / Screw	Multi-Unit Abutment / Screw	Multi-Unit Prosthetic Screw
3.0 mm	15 Ncm	15 Ncm	15 Ncm
3.5 mm, 4.5 mm	35 Ncm	30 Ncm	15 Ncm

NOTE: Any screw-retained prosthetic component not listed in the table above should be hand-tightened only. These values are accurate as of the time printing. Always verify torque values with the implant manufacturer.

■ INSTRUCTIONS FOR USE — INCLUSIVE® ABUTMENTS

Inclusive Abutments are prefabricated, screw-retained intraoral abutments intended to be connected directly to an endosseous implant for retention of a cemented dental prosthesis. They may be indicated for single- and multiple-tooth restorations. Titanium abutments are machined from titanium alloy and attached to the implant fixture with an abutment screw compatible with the restorative instrumentation of the specified implant system.

Capture Implant Placement

Take an implant-level impression utilizing the preferred technique (direct, indirect, or intraoral scan). Submit the impression to the laboratory.

CAD/CAM Preparation

Laboratory — Design the Restoration

- 1) Create a soft tissue model from an implant-level impression.
- 2) Select the appropriate laboratory scanning abutment to capture the implant angulation, position, and abutment connection orientation. Follow manufacturer instructions to obtain all necessary scans to construct an accurate, complete 3-D model.
- 3) Design the abutment according to the patient's clinical needs, taking care to ensure adequate support for the eventual restoration, including appropriate interproximal and occlusal space. Produce a digital design file.
- 4) Send the digital design file to a milling center to manufacture the patient-specific implant abutment.

Milling Center — Fabricate the Restoration

- 1) Select the appropriate Inclusive® Abutment Blank based on the system and platform size of the implant seated in the patient's mouth.
- 2) Fabricate the restoration using CAD/CAM techniques. Veneer as necessary. If a screw-retained hybrid restoration is indicated, fabricate the superstructure (i.e., zirconia coping or crown) and lute it to the titanium abutment. The superstructure is to be bonded to the titanium abutment using MonoCem® Self-Adhesive Resin Cement (Shofu Dental Corporation; San Marcos, Calif.).

Non-CAD/CAM Preparation

Laboratory — Fabricate the Restoration

- 1) Create a soft tissue model from an implant-level impression.
- 2) Select the appropriate Inclusive Abutment based on the system, platform size, location, and occlusal clearance of the implant seated in the patient's mouth.
- 3) Seat the abutment completely into the implant analog on the working model, making sure that the anti-rotational features of the connection interface are fully engaged and the contours of the emergence profile (if applicable) are esthetically oriented.
- 4) Insert the appropriate compatible abutment screw into the abutment's screw access hole and hand-tighten using the appropriate driver.
- 5) Fabricate the restoration using conventional casting techniques. Veneer as necessary. If a screw-retained hybrid restoration is indicated, lute the ceramic crown to the titanium abutment. The ceramic crown is to be bonded to the titanium abutment using MonoCem Self-Adhesive Resin Cement (Shofu Dental Corporation; San Marcos, Calif.).

Manual Adjustment

NOTE: Due to the high thermal conductivity of titanium, titanium abutments should not be modified in the oral cavity. Any necessary modifications should be made extraorally.

- 1) Seat the abutment completely into an implant analog retained by an analog holder or the implant analog captured in the working model, making sure that the anti-rotational features of the connection interface are fully engaged and the contours of the emergence profile (if applicable) are esthetically oriented.
- 2) Insert the appropriate compatible abutment screw into the abutment's screw access hole and hand-tighten using the appropriate driver.
- 3) Using a fine-diamond or carbide bur, modify the abutment as needed.
- 4) With a silicone-based rubber wheel or point, refine the abutment along the margins.

Deliver the Final Restoration

- 1) Seat the titanium abutment or screw-retained hybrid restoration completely into the implant, making sure that the anti-rotational features of the connection interface are fully engaged and the contours of the sculpted emergence profile are esthetically oriented.
- 2) Insert the appropriate compatible abutment screw into the screw access hole and hand-tighten using the appropriate driver. It is strongly recommended that a radiograph of the connection site be taken to confirm complete seating of the abutment or hybrid restoration before proceeding.
- 3) Using the appropriate driver in conjunction with a properly metered torque wrench, tighten the abutment or hybrid restoration to the implant manufacturer's recommended torque value.
- 4) Fill the screw access hole with Teflon tape, gutta-percha, or other suitable material.
- 5) If the restoration is of a screw-retained hybrid design, cover the screw access hole with flowable composite, and cure. Otherwise, follow applicable cementation procedures to affix the definitive restoration to the abutment.

■ INSTRUCTIONS FOR USE — INCLUSIVE® TITANIUM ESTHETIC ABUTMENTS

Inclusive Titanium Esthetic Abutments are prefabricated, screw-retained intraoral abutments intended to be connected directly to an endosseous implant for retention of a cemented dental prosthesis. They may be indicated for single- and multiple-tooth restorations. Titanium esthetic abutments are precisely machined from titanium alloy and attached to the implant fixture with a titanium screw. Unlike the circular emergence profile of standard stock abutments, esthetic abutments are manufactured with a tapered emergence profile for more natural-looking contouring of the soft tissue at the implant site. Each esthetic abutment is specific to the restorative platform of the implant, and anatomically designed for the connection site's region on the ridge (anterior or posterior). In addition to the standard, straight abutment body, angled abutment bodies, produced with a 15 degree slope of one hemisphere to compensate for an undesirable path of insertion resulting from excessive implant angulation, are available.

Capture Implant Placement

Take an implant-level impression utilizing the preferred technique (direct, indirect, or intraoral scan). Submit the impression to the laboratory.

Laboratory — Fabricate the Restoration

- 1) Follow pouring procedures for the appropriate die stone to produce a working model and articulate with a bite registration.
- 2) Select the appropriate Inclusive Titanium Esthetic Abutment based on the system, platform size, location, angulation, and occlusal clearance of the implant seated in the patient's mouth.
- 3) Seat the abutment completely into the implant analog on the working model, making sure that the anti-rotational features of the connection interface are fully engaged and the contours of the sculpted emergence profile are esthetically oriented. For angled abutments, the tapered side should be oriented nearest vertical along the same plane as the implant.
- 4) Insert an Inclusive Titanium Screw into the abutment's screw access hole and hand-tighten using the appropriate driver.
- 5) Fabricate the restoration using conventional casting or CAD/CAM techniques. Veneer as necessary. If a screw-retained hybrid restoration is indicated, lute the ceramic crown to the titanium abutment.

Deliver the Final Restoration

- 1) Seat the titanium esthetic abutment or screw-retained hybrid restoration completely into the implant, making sure that the anti-rotational features of the connection interface are fully engaged and the contours of the sculpted emergence profile are esthetically oriented. For angled abutments, the tapered side should be oriented nearest vertical along the same plane as the implant.
- 2) Insert an Inclusive Titanium Screw into the screw access hole and hand-tighten using the appropriate driver. It is strongly recommended that a radiograph of the connection site be taken to confirm complete seating of the abutment or hybrid restoration before proceeding.
- 3) Using the appropriate driver in conjunction with a properly metered torque wrench, tighten the abutment or hybrid restoration to the implant manufacturer's recommended torque value.
- 4) Fill the screw access hole with Teflon tape, gutta-percha, or other suitable material.
- 5) Follow applicable cementation procedures to affix the definitive restoration to the abutment. Or, if the restoration is of a screw-retained hybrid design, cover the screw access hole with flowable composite, and cure.

■ INSTRUCTIONS FOR USE — INCLUSIVE® MULTI-UNIT ABUTMENTS

Inclusive Multi-Unit Abutments are prefabricated, screw-retained intraoral abutments intended to be connected directly to endosseous implants in partially or fully edentulous patients for the retention of cast or milled bar overdentures. Multi-unit abutments are machined from titanium alloy, and are available with a variety of collar heights to achieve optimal emergence from shallow or deep gingival wells. Each Inclusive Multi-Unit Abutment is delivered sterile, with a carrier color-coded to indicate the restorative platform of the implant.

Straight Multi-Unit Abutments lack any anti-rotational features at the implant-abutment interface. The apical portion of a straight Multi-Unit Abutment is threaded for integration with the internal cavity of the implant. For abutment delivery, the occlusal surface features a male hex head compatible with the Inclusive Multi-Unit Driver.

Angled Multi-Unit Abutments of 17 degrees or 30 degrees enable clinicians to compensate for the divergence of implants or to otherwise accommodate an angled path of insertion. Angled Multi-Unit Abutments feature an anti-rotational connection interface specific to the matching implant platform, and are attached to the implant fixture with an angled Multi-Unit Abutment screw compatible with the restorative instrumentation of the specified implant system. Both straight and angled Multi-Unit Abutments feature a female connection port at the coronal apex, to allow for the attachment of a screw-retained or fixed-removable dental prosthesis with a multi-unit restorative screw (Inclusive Prosthetic Screw).

The axial tilt of an Inclusive Angled Multi-Unit Abutment (angular divergence from path of insertion) is designed and manufactured to lie along a *plane* of the implant connection geometry, as opposed to a corner or junction. To maximize the angle-correcting attributes of the Multi-Unit Abutment, be sure to rotate the implant upon final seating so that one side of the internal connection geometry (flat or lobe) is oriented to serve as the base of angulation, in accordance with the restorative treatment plan.

Place the Multi-Unit Abutment

- 1) Select the appropriate Inclusive Multi-Unit Abutment based on platform size, implant angle, and depth of the soft-tissue well.
- 2) Retrieve the abutment from its packaging. To maintain the sterility of the Multi-Unit Abutment, be careful to handle only by the carrier.
- 3) (a) *For Straight Abutments:* Using the carrier, seat the abutment into the implant and hand-tighten. Remove the carrier by pulling the apex of the carrier toward the facial. (b) *For Angled Abutments:* Using the carrier, seat the abutment into the implant until the anti-rotational features of the connection interface are engaged. Lift and rotate as necessary to orient the angle in the required direction. Hand-tighten the Inclusive Angled Multi-Unit Abutment Screw using the appropriate driver. Twist the carrier counterclockwise to unscrew the carrier from the abutment.

NOTE: It is strongly recommended that a radiograph of the connection site be taken to confirm complete seating of the abutment before proceeding.

- 4) Using the appropriate driver in conjunction with a properly metered torque wrench, tighten the Multi-Unit Abutment or Angled Multi-Unit Abutment Screw to the implant manufacturer's recommended torque value.

Passive Temporization of Multi-Unit Abutments

- 1) If delayed loading is desired, cover each Inclusive Multi-Unit Abutment with an Inclusive Multi-Unit Temporary Healing Cap and hand-tighten with the Inclusive Prosthetic Screw provided, using the appropriate driver. Do not overtighten.
- 2) Using the patient's existing denture or other prosthesis, relieve the area directly above the placement of each temporary healing cap until the denture rests on the ridge.
- 3) Follow procedures to reline the denture over the temporary healing caps, using tissue conditioner or soft relining material. The temporized denture can be used during healing.

NOTE: For a temporization technique involving loading, please refer to the *Inclusive Prosthetic Components Restorative Manual*.

Capture Multi-Unit Abutment Placement

Take an abutment-level impression utilizing the preferred technique (direct, indirect, or intraoral scan). Submit the impression to the laboratory for the fabrication of a working cast and verification index.

Denture Protocol

Follow appropriate denture protocol in accordance with the patient-specific treatment plan. When trying in the various setups (e.g., verification index, occlusal rim, wax setup, retention bar), hand-tighten to the Multi-Unit Abutments with prosthetic screws, using the appropriate driver. Start from the distal and move forward, alternating between sides of the ridge. Always confirm complete, passive seating, modifying the setup as needed.

Deliver the Final Restoration

- 1) Remove any temporary prosthesis.
- 2) Confirm that each Multi-Unit Abutment is tightened to the implant manufacturer's recommended torque value.
- 3) Align the prosthesis onto the abutments. Beginning with the midmost screw access channel, hand-tighten an Inclusive Prosthetic Screw into the Multi-Unit Abutment. Repeat for each abutment, working outward and alternating left to right.
- 4) Confirm appropriate seating. With the same middle-out, left-to-right technique, tighten each prosthetic screw to the recommended torque value.
- 5) Check comfort and occlusion, and make any necessary adjustments.
- 6) Place Teflon tape, gutta-percha, or other suitable material.
- 7) Cover each screw access hole with flowable composite and cure.

SURGICAL INSTRUMENTS

■ Description

Inclusive surgical instruments and surgical/restorative accessories are made out of the following materials: titanium alloy, gold alloy, polymers, and stainless steel. They are designed for use with Inclusive dental implants and prosthetic components. For specific product identification and contents, please refer to individual component packaging and appropriate product catalogs and/or user manuals.

■ Processing/Reprocessing Instructions

Autoclavable Cassette and surgical instruments are delivered non-sterile and must be cleaned and sterilized before the first use, after every use, and every time they are contaminated. Use of non-sterile devices may lead to infection of tissue or infectious diseases. If sterile barrier has been compromised, the device must be re-sterilized prior to use. Automated cleaning may not be effective. The manual cleaning process below is validated and recommended.

Point of Use:

- 1) After surgery, collect used instruments separately, and place them in a sterile water bath to prevent contamination from drying. Do not place soiled instruments in the tray.
- 2) Use an enzymatic spray or wipe to prevent contamination from drying on the Autoclavable Cassette.

Containment and Transportation:

Follow the clinical internal procedures for the transportation of contaminated devices. Perform cleaning as soon as possible. When longer delays are expected, immerse the devices in a sterile water bath to avoid drying of debris. Delay in reprocessing must be kept to a minimum to avoid contaminants drying. Prolonged storage in a sterile water bath may result in degradation of the product.

Pre-Treatment:

- 1) Remove all instruments from the Autoclavable Cassette. Properly dispose of any damaged instruments.
- 2) Disassemble the Autoclavable Cassette to Lid, Tray, and Base.

Manual Cleaning for Autoclavable Cassette:

1. Rinse the lid, tray, and base of the Autoclavable Cassette under running, cold utility water for a minimum of 1 minute to remove excess contamination.
2. Prepare an enzymatic detergent per the manufacturer's recommendations using utility water.
3. Fully immerse all components of the Autoclavable Cassette in the prepared detergent bath and allow them to soak for minimum duration recommended by the manufacturer.
4. After the soak, but while still immersed, use a clean, soft-bristled brush to brush all components of the Autoclavable Cassette to remove visible soil. Use a soft bristled lumen/tube brush to remove visible soil from the holes of the Autoclavable Cassette, as needed.
 - Ensure to brush all hard-to-reach areas, such as the parts between the grommets. Remove the grommet if visible soil is present between the grommet and the tray interface by pushing the grommets up from the bottom. For offset grommets, push down from the top.
5. Remove all components of the Autoclavable Cassette from the prepared detergent bath and rinse them under running Critical water for a minimum of 1 minute to remove detergent residuals.
6. Dry the Autoclavable Cassette and components using lint-free cloths and/or filtered pressurized air.
7. Inspect the Autoclavable Cassette and components for visible contamination. Repeat steps 1 through 6 if contamination is observed.

Manual Cleaning for Instruments:

1. If present, special Instructions for Adjustable Torque Wrench: Disassemble and clean the Adjustable Torque Wrench as per manufacturer's instructions for use.
2. Rinse gross debris for a minimum of thirty (30) seconds under running cold utility water.
3. While completely submerged, brush the exterior of device with a soft bristled brush to remove visible debris for a minimum of 20 seconds. If applicable, clean the interior lumina of instrument with lumen/pipe brushes.
 - Flush the inner surfaces, lumina, and cavities (where applicable) with 20 ml cold utility water using an irrigation needle connected to a disposable syringe.
 - While completely submerged, brush the inner surfaces, lumina and cavities (where applicable) with appropriately sized lumen/pipe brush for a minimum of 20 seconds until all visible soil is removed.
4. Inspect instruments visually and check for cleanliness, and for visible damage and/or wear. Discard damaged or worn instrument. If any remaining bone fragments, visible soil, or residual debris remain repeat the scrubbing and rinsing procedure as necessary.
5. Prepare an enzymatic solution for ultrasonic cleaner per the manufacturer's recommendations.
6. Place instruments in a single layer into the ultrasonic bath. Make sure the instruments are completely submerged and do not touch one another. Turn on ultrasonic bath for five (5) minutes. Note: process Titanium instruments separately from Stainless Steel instruments.

7. Remove the device and thoroughly rinse all surfaces of device with tepid running critical water for a minimum of ten (10) seconds to remove all cleaning agent.
8. Flush the inner surfaces, lumina, and cavities (where applicable) with 20 ml of tepid critical water using a new irrigation needle connected to a new disposable syringe.
9. Dry instrument using lint-free cloths and/or filtered pressurized air.
10. If the instrument is still dirty, repeat the cleaning steps 1 through 9 above.

The recommended sterilization process is based on the ANSI/AAMI/ISO 17665-1 and ANSI/AAMI ST79 guidelines, as follows:

- **Sterilization:** Gravity-fed sterilizers: Autoclave in sterilization pouch for 15 minutes at 132°C (270°F). Allow sterilized components to dry for at least 30 minutes.

NOTE: The validated procedures require the use of FDA-cleared sterilization trays, wraps, biological indicators, chemical indicators, and other sterilization accessories labeled for the sterilization cycle recommended. The healthcare facility should monitor the sterilizer for the facility according to an FDA-recognized sterility assurance standard such as ANSI/AAMI ST79.

■ Warnings

Prior to surgery, ensure that instruments and accessories are complete, functional, and available in the correct quantities. Instruments may be used for up to five preparations. For best results, replace regularly.

■ Precautions

For best results, please observe the following precautions:

- Proper surgical protocol should be strictly adhered to.
- All surgical instruments used must be in good condition and should be used carefully to avoid damage to implants or other components.
- Since implant components and their instruments are very small, precautions should be taken to ensure that they are not swallowed or aspirated by the patient.

SYMBOLS GLOSSARY

Symbol	Symbol Ref. No.	Symbol Title	Designation No.	Explanatory Text
	5.2.4	Sterile with Gamma Radiation	EN ISO 15223-1	This symbol indicates that this device has been sterilized using irradiation.
	5.2.8	Do Not Use if Package is Damaged	EN ISO 15223-1	This symbol indicates that this device should not be used if the package has been damaged or opened.
	5.2.7	Non-Sterile	EN ISO 15223-1	This device has not been subjected to a sterilization process.
	5.4.2	Do not Re-use	EN ISO 15223-1	This device is intended for one use, or for use on a single patient during a single procedure.
	5.2.6	Do not Resterilize	EN ISO 15223-1	This symbol indicates that this device is not to be resterilized.
	5.1.4	Use-by Date	EN ISO 15223-1	This symbol indicates the date (YYYY-MM-DD) after which this device is not to be used.
	Sec. 801.109(b)(1)	By Prescription Only	21 CFR Part 801	Caution: Federal law restricts this device to sale by, or on the order of, a licensed dentist or physician.
	5.1.6	Catalog Number	EN ISO 15223-1	This symbol indicates PrismaTik Dentalcraft's catalog number so that this device can be identified.
	5.1.5	Lot/Batch Number	EN ISO 15223-1	This symbol indicates PrismaTik Dentalcraft's lot/batch number so that the lot/batch of this device can be identified.
	5.4.3	Consult Instructions For Use	EN ISO 15223-1	This symbol indicates the need of the user to consult the instructions for use.
	5.1.1	Manufacturer Date of Manufacture (YYYY-MM-DD)	EN ISO 15223-1	This symbol indicates the manufacturer and the date of manufacture of this device.

Made in USA

Within the U.S.: 888-303-3975

Within Canada: 888-278-0414

glidewell.com

PrismaTik Dentalcraft, Inc.
2144 Michelson Drive
Irvine, CA 92612, USA

Instructions for Use

IMPORTANT INFORMATION — PLEASE READ

Caution: U.S. federal law restricts this device to sale by, or on the order of, a licensed dentist or physician.

■ General Information

The Inclusive Mini Implant System consists of dental implants, prosthetic components, surgical instrumentation, and related accessories packaged under the Inclusive brand name for use by qualified, licensed clinicians and laboratory technicians fully trained in their application.

For specific product identification and contents, please refer to individual product labels and the following catalog:

- Inclusive Dental Implant System Product Catalog (**MKT 787**)

For detailed information on the specifications and intended use of a particular product, please refer to the following user manual:

- Inclusive Dental Implant System Surgical Manual (**UM 4236**)

■ Explanation of Label Symbols

The symbols glossary is provided on page 15 of this IFU document.

■ Disclaimer of Liability

The guidelines presented herein are not adequate to allow inexperienced clinicians to administer professional implant treatment or prosthetic dentistry, and are not intended to substitute for formal clinical or laboratory training. Inclusive devices should only be used by individuals with training and experience specific to their clinically accepted application.

Prismatik Dentalcraft, Inc. is not liable for damages resulting from treatment outside of our control. The responsibility rests with the provider.

■ MRI

The Inclusive Mini Implant System has not been evaluated for safety and compatibility in the magnetic resonance (MR) environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the Inclusive Mini Implant System in the MR environment is unknown. Scanning a patient who has the device may result in patient injury.

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■ Description

Inclusive Mini Implants are endosseous devices manufactured from titanium alloy. They are compatible with the prosthetic components and surgical instrumentation of the Inclusive Mini Implant System.

■ Indications for Use

Inclusive Mini Implants are self-tapping threaded titanium screws indicated for long-term applications. Inclusive Mini Implants may also be used for provisional applications. These devices will allow immediate loading and long-term stabilization of dentures and provisional stabilization of dentures while standard implants heal. To be used for immediate loading only in the presence of primary stability and appropriate occlusal loading.

■ Contraindications

Patients should be evaluated before the time of surgery for factors that put them at risk from the implant placement procedure, or that may affect healing of bone or surrounding soft tissue.

- Implant placement in patients medically unfit for oral surgical procedures is contraindicated.
- Patients with systemic, localized or pharmaceutical treatment factors that compromise their ability to heal should be carefully evaluated.
- Do not place Inclusive Mini Implants if there is not adequate bone width or height to contain the implant.
- Do not reuse Inclusive Mini Implants. Reuse of Inclusive Mini Implants are not indicated for abutment or crown restorations.

■ Warnings

- Inclusive Mini Implants are intended to be used on an individual patient only. The reuse of such device on another patient is not recommended due to the risks of cross-contamination or infection.
- Inclusive Mini Implants may only be used for their intended purpose in accordance with general rules for dental/surgical treatment, occupational safety, and accident prevention. They must only be used for dental procedures with the restorative components they were designed for. If the indications and intended use are not clearly specified, treatment should be suspended until these considerations have been clarified.
- The following instructions are not sufficient to allow inexperienced clinicians to administer professional prosthetic dentistry. Inclusive Mini Implants, surgical instruments, and restorative components must only be used by dentists and surgeons with training/experience with oral surgery, prosthetics and biochemical requirements, as well as diagnosis and preoperative planning.
- The implant site should be inspected for adequate bone by radiographs, palpations and visual examination. Determine the location of nerves and other vital structures and their proximity to the implant site before any drilling to avoid potential injury, such as permanent numbness to the lower lip and chin.
- Prior to surgery, ensure that the needed components, instruments and help materials are complete, functional and available in the correct amounts.
- Treatment of children is not recommended until growth is finished and epiphyseal closure has occurred.
- Inclusive Mini Implants should always be used in sufficient quantity to prevent excessive stress on the implants; at least one pair in all cases. Absolute success cannot be guaranteed. Factors such as infection, disease, and inadequate bone quality and/or quantity can result in osseointegration failures following surgery or initial osseointegration.
- Inclusive Mini Implants can distort images obtained via magnetic resonance imaging (MRI).
- Prismatik Dentalcraft, Inc. is not liable for damages resulting from treatment outside of our control. The responsibility rests with the provider.

■ Precautions

Surgical Procedures

Minimizing tissue damage is crucial to successful implant osseointegration. In particular, care should be taken to eliminate sources of infection, contaminants, surgical and thermal trauma. Risk of osseointegration failure increases as tissue trauma increases. For best results, please observe the following precautions:

- All drilling procedures should be performed at 2000 RPM or less under continual, copious irrigation.
- All surgical instruments used must be in good condition and should be used carefully to avoid damage to implants or other components. Implants should be placed with sufficient stability; however, excessive insertion torque may result in implant fracture, or fracture or necrosis of the implant site. The proper surgical protocol should be strictly adhered to.
- Since implant components and their instruments are very small, precautions should be taken to ensure that they are not swallowed or aspirated by the patient.
- Prior to surgery, ensure that the needed components, instruments and ancillary materials are complete, functional and available in the correct quantities.

Prosthetic Procedures

Following successful placement of Inclusive Mini Implants, verify primary stability and appropriate occlusal loading before proceeding with the placement of a permanent or provisional prosthesis. All components that are used intraorally should be secured to prevent aspiration or swallowing. Distribution of stress is an important consideration. Care should be taken to avoid excessive loads significantly transverse to the implant axes.

■ Sterility

Inclusive Mini Implants are shipped sterile. They should not be resterilized. They are for single use only, prior to the expiration date. Do not use implants if the packaging has been compromised or previously opened.

■ Storage and Handling

Inclusive Mini Implants must be stored in a dry location at room temperature, in their original packaging. Inclusive Mini Implants are packaged sterile. Do not handle implant surfaces directly. Users are advised to visually inspect packaging to ensure seals and contents are intact prior to use. Please refer to the individual product label for all relevant product information and cautions.

■ INSTRUCTIONS FOR USE

Case Planning

Prescribed implant length and diameter should take into account crestal width, cortical thickness, bone density, and any other relevant clinical factors. Use appropriate radiography in mandibular cases to identify the location of the inferior alveolar nerve, including a possible anterior loop. Use appropriate radiography in maxillary cases to identify the location of the sinuses. When patient evaluation is complete, establish the number of Inclusive Mini Implants required for denture stabilization and identify the appropriate implant sites. In mandibular cases, it is recommended that four (4) mini implants be placed within the symphysis area with as wide an anterior-posterior spread as possible while still ensuring an adequate margin of safety from the nerve. In maxillary cases, it is recommended that six (6) mini implants be placed anterior to the sinuses. Wider implants are often preferred for softer bone types. Mini implants should be placed with a minimum of 7 mm between implants, to accommodate the size of the O-ring housings. The housings can accommodate up to 30 degrees of angular divergence between mini implants. However, implants should be placed as parallel to one another as possible to provide ideal prosthetic fit and to avoid excessive wearing of the O-rings. There should be at least 8 mm of vertical space from the top of the implant collar to allow for adequate thickness of the prosthesis. The denture teeth would be in addition to this space.

Drilling Protocol

Mark each implant site on the patient's tissue. Select the appropriate Cortical Bone Drill (1.5 mm, 1.7 mm, or 2.4 mm), as determined by the patient's bone density and the diameter of the implant to be placed. Carefully place the drill directly above the implant site and gently drill through the tissue and alveolar crest using an in-and-out motion and profuse, sterile irrigation to a depth of one-third (1/3) to one-half (1/2) the length of the implant threads. If placing 3.0 mm diameter Inclusive Mini Implants, continue drilling to a depth of at least two-thirds (2/3) the length of the implant threads. For the majority of implant sites, this is the extent of the drilling that is required. However, in dense bone, the drilling depth may need to be greater. The goal is to achieve high primary stability with an insertion torque of approximately 35 Ncm, taking care not to exceed the recommended maximum of 45 Ncm.

Implant Placement

Remove the titanium implant holder from its packaging and place it onto a sterile field. Use slight finger pressure to pinch the occlusal end of the implant in its holder while inserting the appropriate Mini Implant Driver. Gently rotate implant and holder, allowing the driver to engage the implant connection. With the driver securely attached to the implant, squeeze the opposing end of the holder to disengage the implant from the holder. Transport the implant to the prepared site, and insert into the osteotomy. Rotate clockwise with applied pressure to engage the self-tapping threads. Avoid lateral forces, which can affect the angulation and final alignment of the implant.

NOTE: Apply pressure to ensure the driver is fully engaged with the implant prior to disengaging the titanium holder.

Final Insertion

With the implant threaded securely in its proper site, slide the Ratchet Wrench fully into place over the mini implant driver. Turn the wrench clockwise in small increments of approximately 90 degrees, pausing between rotations to allow the bone to expand. Avoid lateral forces, which can affect the final angulation of the implant. Optimal final insertion of the implant leaves the implant head fully exposed, while the collar is embedded in the gingiva with no threads visible. For immediate loading of the implant, final torque at seating should be 30–35 Ncm minimum. Exceeding 45 Ncm torque during implant placement is not recommended.

NOTE: If the implant cannot be fully seated using the recommended torque, it may be necessary to reverse the implant from the site and drill again to increase the depth of the osteotomy. For positive long-term prognosis, solid resistance must be met during final insertion. Inadequate resistance contraindicates primary stability and loading. In such instances, a larger implant should be placed, or a new implant site determined.

Impression Procedure

An impression procedure is required whenever a new removable prosthesis is going to be fabricated. Based on the clinician's preference, the O-ring housings can be processed into the denture, or space made and the housings picked up chairside.

Step 1: Seat the Copings – Snap a Mini Implant Impression Coping onto the head of each Inclusive Mini Implant. If gingival tissue prevents full engagement of a coping onto an implant, take an impression of the mini implant without the use of impression copings, or trim the tissue.

Step 2: Seat the Impression – Standard impression techniques are used to pick up the impression copings, recording each implant’s position easily and accurately.

Step 3: Remove the Impression – Once the impression has fully set, carefully remove the tray from the patient’s mouth and verify that all impression copings have been captured accurately in the impression.

Step 4: Insert the Replicas – This step can be performed in the clinic or at the dental laboratory. Align the squared neck of a Mini Implant Replica with the squared opening at the base of the impression coping. Press the replica into the coping until it snaps into proper position. Insert a replica into each coping and prepare the impression to be used to fabricate a stone model.

Step 5: Fabricate the Model – Use standard laboratory procedure to fabricate a soft tissue model.

Soft Denture Reline

A soft denture reline procedure is used when immediate loading with the O-rings is contraindicated, as in the case of a transitional prosthesis, or whenever the Inclusive Mini Implants are placed in soft bone (such as the maxilla or a Type III mandible). Following an appropriate healing period, the soft inner liner can be replaced with a hard pick-up of the O-ring housings to increase the level of retention.

Step 1: Prepare the Denture – Relieve the patient’s existing denture to make room for the implant heads. The positions of the implants can be identified using a color transfer applicator, or by lining the intaglio surface of the denture with impression or bite registration material. An acrylic bur can then be used to relieve the denture. The denture must be sufficiently relieved to seat passively, without resting on or against the implant heads. Lightly roughen the tissue-facing surface of the denture with an acrylic bur, and degrease the surface with isopropyl alcohol.

Step 2: Line the Denture – Apply the selected soft reline material onto the tissue-facing surface of the denture. Seat the denture in the patient’s mouth. Instruct the patient to close with normal pressure into centric occlusion. Allow the soft reline material to set.

Step 3: Final Preparation – Remove the denture from the patient’s mouth and trim excess material with fine scissors or a surgical blade. Do not remove the palate of a maxillary denture during this stage. Instruct the patient to keep the denture in place for the first 48 hours following placement, to prevent gingival overgrowth.

Hard Denture Reline

A hard denture reline procedure is used to incorporate the retention caps (O-ring Housings) that cover the Inclusive Mini Implants in the patient’s final prosthesis. This loading procedure can typically be performed immediately after placement of the Inclusive Mini Implants, provided primary stability and appropriate occlusal loading are assured. Primary stability is generally indicated when 35 Ncm of torque resistance is achieved, with implants seated at the appropriate gingival depth.

Step 1: Prepare the Denture – Mark the location of the implants on the intaglio surface of the patient’s existing denture. This can be done using a color transfer applicator, or by lining the intaglio surface of the denture with impression or bite registration material. Relieve the denture to make room for the O-ring housings. This can be done by creating a space for each housing where marked (or by burring a full trough). The denture must be sufficiently relieved to seat passively, without resting on or against the implant heads.

Step 2: Block Out the Implant Heads – Use a rubber dam or trim the Blockout Shims to the appropriate length in order to completely mask the exposed neck of each implant beneath the O-ball head. This is critical to prevent pick-up material from flowing under the O-ball. Place an O-ring housing on each mini implant, checking for passive fit over the blockout shims. Place the denture in the patient’s mouth, checking for passive fit over implants and housings.

Step 3: Line the Denture – Apply a thin layer of adhesive on the intaglio surface of the denture. Place hard pick-up material directly onto the O-ring housings and into the housing spaces (or trough) in the denture. Seat the denture in the patient’s mouth. Instruct the patient to close with normal pressure into centric occlusion. Allow the hard pick-up material to set.

Step 4: Final Preparation – Remove the denture and all blockout shims. Trim and polish. Instruct the patient to keep the denture in place for the first 48 hours following implant placement, to prevent gingival overgrowth.

PROSTHETIC COMPONENTS

■ Description

Prosthetic components of the Inclusive Mini Implant System, consisting of O-ring housings, analogs, impression copings, and related restorative accessories, are manufactured from titanium alloy or polymers. They are compatible with Inclusive Mini Implants. For product-specific descriptions and sterility information, please refer to the individual product labels and appropriate Inclusive product catalog and/or user manual.

■ Warnings

Inclusive Mini Implant prosthetic components are intended to be used on an individual patient only. The reuse of such device on another patient is not recommended due to the risks of cross-contamination or infection.

■ Adverse Effects

The following adverse effects have been observed when using prosthetic components and accessories:

- Components used in the patient’s mouth have been aspirated or swallowed.
- Components are not fully seated due to inadequate application of force.

■ Precautions

Inclusive Mini Implant prosthetic components may only be used for their intended purpose in accordance with general rules for dental/prosthetic treatment, occupational safety, and accident prevention. They must only be used for dental procedures with the Inclusive Mini Implant System. If the indications and intended use are not clearly specified, treatment should be suspended until these considerations have been clarified. All components that are used intraorally must be secured to prevent aspiration or swallowing. Prior to placement, ensure that the required components, instruments, and ancillary materials are complete, functional, and available in the correct quantities.

■ Side Effects

No side effects, according to current knowledge.

■ Sterility

Inclusive Mini Implant prosthetic components are shipped non-sterile. Non-sterile components must be cleaned, disinfected, and sterilized prior to clinical use, according to a validated method.

- **Cleaning:** Wash using a broad spectrum cleaning solution, followed by thorough rinsing.

The recommended disinfection process is based on ANSI/AAMI ST79 guidelines, as follows:

- **Disinfection:** Immerse abutments in disinfectant¹, rinse with distilled water.

The recommended sterilization process is based on the ANSI/AAMI/ISO 17665-1 and ANSI/AAMI ST79 guidelines, as follows:

- **Sterilization:** Gravity-fed sterilizers: Autoclave in sterilization pouch for 30 minutes at 121°C (250°F). Devices are to be used immediately after sterilization.

NOTE: The validated procedures require the use of FDA-cleared sterilization trays, wraps, biological indicators, chemical indicators, and other sterilization accessories labeled for the sterilization cycle recommended. The healthcare facility should monitor the sterilizer for the facility according to an FDA-recognized sterility assurance standard such as ANSI/AAMI ST79.

¹Oral disinfectant containing *Chlorhexidine* is recommended; refer to the disinfectant manufacturer's instructions.

SURGICAL INSTRUMENTS

■ Description

Inclusive surgical instruments and surgical/restorative accessories are made out of the following materials: titanium alloy, gold alloy, polymers, and stainless steel. They are designed for use with Inclusive dental implants and prosthetic components. For specific product identification and contents, please refer to individual component packaging and appropriate product catalogs and/or user manuals.

■ Processing/Reprocessing Instructions

Autoclavable Cassette and surgical instruments are delivered non-sterile and must be cleaned and sterilized before the first use, after every use, and every time they are contaminated. Use of non-sterile devices may lead to infection of tissue or infectious diseases. If sterile barrier has been compromised, the device must be re-sterilized prior to use. Automated cleaning may not be effective. The manual cleaning process below is validated and recommended.

Point of Use:

- 1) After surgery, collect used instruments separately, and place them in a sterile water bath to prevent contamination from drying. Do not place soiled instruments in the tray.
- 2) Use an enzymatic spray or wipe to prevent contamination from drying on the Autoclavable Cassette.

Containment and Transportation:

Follow the clinical internal procedures for the transportation of contaminated devices. Perform cleaning as soon as possible. When longer delays are expected, immerse the devices in a sterile water bath to avoid drying of debris. Delay in reprocessing must be kept to a minimum to avoid contaminants drying. Prolonged storage in a sterile water bath may result in degradation of the product.

Pre-Treatment:

- 1) Remove all instruments from the Autoclavable Cassette. Properly dispose of any damaged instruments.
- 2) Disassemble the Autoclavable Cassette to Lid, Tray, and Base.

Manual Cleaning for Autoclavable Cassette:

1. Rinse the lid, tray, and base of the Autoclavable Cassette under running, cold utility water for a minimum of 1 minute to remove excess contamination.
2. Prepare an enzymatic detergent per the manufacturer's recommendations using utility water.
3. Fully immerse all components of the Autoclavable Cassette in the prepared detergent bath and allow them to soak for minimum duration recommended by the manufacturer.
4. After the soak, but while still immersed, use a clean, soft-bristled brush to brush all components of the Autoclavable Cassette to remove visible soil. Use a soft bristled lumen/tube brush to remove visible soil from the holes of the Autoclavable Cassette, as needed.
 - Ensure to brush all hard-to-reach areas, such as the parts between the grommets. Remove the grommet if visible soil is present between the grommet and the tray interface by pushing the grommets up from the bottom. For offset grommets, push down from the top.
5. Remove all components of the Autoclavable Cassette from the prepared detergent bath and rinse them under running Critical water for a minimum of 1 minute to remove detergent residuals.
6. Dry the Autoclavable Cassette and components using lint-free cloths and/or filtered pressurized air.
7. Inspect the Autoclavable Cassette and components for visible contamination. Repeat steps 1 through 6 if contamination is observed.

Manual Cleaning for Instruments:

1. If present, special Instructions for Adjustable Torque Wrench: Disassemble and clean the Adjustable Torque Wrench as per manufacturer's instructions for use.
2. Rinse gross debris for a minimum of thirty (30) seconds under running cold utility water.
3. While completely submerged, brush the exterior of device with a soft bristled brush to remove visible debris for a minimum of 20 seconds. If applicable, clean the interior lumina of instrument with lumen/pipe brushes.
 - Flush the inner surfaces, lumina, and cavities (where applicable) with 20 ml cold utility water using an irrigation needle connected to a disposable syringe.
 - While completely submerged, brush the inner surfaces, lumina and cavities (where applicable) with appropriately sized lumen/pipe brush for a minimum of 20 seconds until all visible soil is removed.
4. Inspect instruments visually and check for cleanliness, and for visible damage and/or wear. Discard damaged or worn instrument. If any remaining bone fragments, visible soil, or residual debris remain repeat the scrubbing and rinsing procedure as necessary.
5. Prepare an enzymatic solution for ultrasonic cleaner per the manufacturer's recommendations.
6. Place instruments in a single layer into the ultrasonic bath. Make sure the instruments are completely submerged and do not touch one another. Turn on ultrasonic bath for five (5) minutes. Note: process Titanium instruments separately from Stainless Steel instruments.
7. Remove the device and thoroughly rinse all surfaces of device with tepid running critical water for a minimum of ten (10) seconds to remove all cleaning agent.
8. Flush the inner surfaces, lumina, and cavities (where applicable) with 20 ml of tepid critical water using a new irrigation needle connected to a new disposable syringe.
9. Dry instrument using lint-free cloths and/or filtered pressurized air.
10. If the instrument is still dirty, repeat the cleaning steps 1 through 9 above.

The recommended sterilization process is based on the ANSI/AAMI/ISO 17665-1 and ANSI/AAMI ST79 guidelines, as follows:

- **Sterilization:** Gravity-fed sterilizers: Autoclave in sterilization pouch for 20 minutes at 132°C (270°F). Devices are to be used immediately after sterilization.

NOTE: The validated procedures require the use of FDA-cleared sterilization trays, wraps, biological indicators, chemical indicators, and other sterilization accessories labeled for the sterilization cycle recommended. The healthcare facility should monitor the sterilizer for the facility according to an FDA-recognized sterility assurance standard such as ANSI/AAMI ST79.

■ Warnings

Prior to surgery, ensure that instruments and accessories are complete, functional, and available in the correct quantities. Instruments may be used for up to five preparations. For best results, replace regularly.

■ Precautions

For best results, please observe the following precautions:

- Proper surgical protocol should be strictly adhered to.
- All surgical instruments used must be in good condition and should be used carefully to avoid damage to implants or other components.
- Since implant components and their instruments are very small, precautions should be taken to ensure that they are not swallowed or aspirated by the patient.

SYMBOLS GLOSSARY

Symbol	Symbol Ref. No.	Symbol Title	Designation No.	Explanatory Text
	5.2.4	Sterile with Gamma Radiation	EN ISO 15223-1	This symbol indicates that this device has been sterilized using irradiation.
	5.2.8	Do Not Use if Package is Damaged	EN ISO 15223-1	This symbol indicates that this device should not be used if the package has been damaged or opened.
	5.2.7	Non-Sterile	EN ISO 15223-1	This device has not been subjected to a sterilization process.
	5.4.2	Do not Re-use	EN ISO 15223-1	This device is intended for one use, or for use on a single patient during a single procedure.
	5.2.6	Do not Resterilize	EN ISO 15223-1	This symbol indicates that this device is not to be resterilized.
	5.1.4	Use-by Date	EN ISO 15223-1	This symbol indicates the date (YYYY-MM-DD) after which this device is not to be used.
	Sec. 801.109(b)(1)	By Prescription Only	21 CFR Part 801	Caution: Federal law restricts this device to sale by, or on the order of, a licensed dentist or physician.
	5.1.6	Catalog Number	EN ISO 15223-1	This symbol indicates PrismaTik Dentalcraft's catalog number so that this device can be identified.
	5.1.5	Lot/Batch Number	EN ISO 15223-1	This symbol indicates PrismaTik Dentalcraft's lot/batch number so that the lot/batch of this device can be identified.
	5.4.3	Consult Instructions For Use	EN ISO 15223-1	This symbol indicates the need of the user to consult the instructions for use.
	5.1.1	Manufacturer Date of Manufacture (YYYY-MM-DD)	EN ISO 15223-1	This symbol indicates the manufacturer and the date of manufacture of this device.

Made in USA

Within the U.S.: 888-303-3975

Within Canada: 888-278-0414

glidewell.com

PrismaTik Dentalcraft, Inc.
2144 Michelson Drive
Irvine, CA 92612, USA

Instructions for Use

IMPORTANT INFORMATION — PLEASE READ

Caution: U.S. federal law restricts this device to sale by, or on the order of, a licensed dentist or physician.

■ Disclaimer of Liability

The guidelines presented herein are not adequate to allow inexperienced clinicians to administer professional implant treatment or prosthetic dentistry, and are not intended to substitute for formal clinical or laboratory training. Inclusive devices should only be used by individuals with training and experience specific to their clinically accepted application.

Prismatik Dentalcraft, Inc. is not liable for damages resulting from treatment outside of our control. The responsibility rests with the provider.

■ Description

Inclusive Prosthetic Components, consisting of abutments, screws, analogs, copings, and related restorative accessories, are manufactured from titanium alloy, gold alloy, polymers, or zirconia. They are for use by qualified, licensed clinicians and laboratory technicians fully trained in their application. Inclusive Prosthetic Components are shipped non-sterile (except for healing abutments and Multi-Unit Abutments).

For product-specific descriptions and sterility information, refer to individual product labels and the Inclusive Prosthetic Components Product Catalog (**MKT 1419**).

For detailed information on the specifications and intended use of a particular product, refer to the Inclusive Prosthetic Components Restorative Manual (**UM 5088**).

■ Indications for Use

Inclusive Abutments are premanufactured prosthetic components directly connected to endosseous dental implants and are intended for use as an aid in prosthetic rehabilitation.

All digitally designed abutments for use with Inclusive Abutments for CAD/CAM are intended to be manufactured at a validated milling center.

Inclusive Multi-Unit Abutments are intended to provide support and retention for multi-unit screw-retained restorations. The 30-degree Multi-Unit Abutments must be used within 45 degrees of parallelism for a splinted restoration. The 17-degree Multi-Unit Abutments must be used within 32 degrees of parallelism for a splinted restoration.

■ Contraindications

Inclusive Abutments

- Wall thickness less than 0.5 mm
- Gingival margin diameter less than 0.5 mm wider than implant
- Gingival height greater than 6 mm**
- Angle corrections of more than 30 degrees
- Angle corrections of more than 20 degrees*
- Less than 0.5 mm margin height
- Less than 4.0 mm abutment post height above the gingival collar

*ONLY applicable to Dentsply Implants Astra Tech Implant System® EV, Dentsply Implants Astra Tech OsseoSpeed™ and Camlog® Screw-Line.

**ONLY applicable to Dentsply Implants Astra Tech Implant System® EV.

Inclusive Hybrid Abutments (titanium bases with zirconia copings)

- Wall thickness less than 0.5 mm
- Gingival margin diameter less than 0.5 mm wider than implant
- Gingival height greater than 6.0 mm
- Angle corrections
- Less than 0.5 mm margin height
- Less than 4.0 mm abutment post height above the gingival collar

NOTE: A patient-specific finished device will consist of both the titanium base and zirconia coping. The titanium base is not intended to provide angle correction. The zirconia coping is intended to be milled straight only.

Inclusive Titanium Esthetic Abutments

- Wall thickness less than 0.5 mm
- Gingival margin diameter less than 0.5 mm wider than implant
- Angle corrections of more than 15 degrees
- Less than 0.5 mm margin height
- Less than 4.0 mm abutment post height above the gingival collar

NOTE: Small-diameter implants with angled abutments are not recommended for the posterior region of the mouth.

Inclusive Multi-Unit Abutments

- Greater than 45 degrees divergence from parallel for a splinted restoration when using 30-degree Multi-Unit Abutments
- Greater than 32 degrees divergence from parallel for a splinted restoration when using 17-degree Multi-Unit Abutments

■ Warnings

An Inclusive abutment is intended to be used on an individual patient only. The reuse of such device on another patient is not recommended due to the risks of cross-contamination or infection. Small-diameter implants with angled abutments are not recommended for the posterior region of the mouth.

■ Adverse Effects

The following adverse effects have been observed when using prosthetic components and accessories:

- Components used in the patient's mouth have been aspirated or swallowed.
- The abutment screw has fractured due to application of excessive torque.

- The abutment is not adequately secured due to inadequate application of torque.

■ Precautions

Inclusive abutments may only be used for their intended purpose in accordance with general rules for dental/prosthetic treatment, occupational safety, and accident prevention. Inclusive abutments must only be used for dental procedures with the implant systems they were designed for. If the indications and intended use are not clearly specified, treatment should be suspended until these considerations have been clarified. All components that are used intraorally must be secured to prevent aspiration or swallowing. Prior to placement, ensure that the required components, instruments, and ancillary materials are complete, functional, and available in the correct quantities.

■ Side Effects

No side effects, according to current knowledge.

■ Storage and Handling

Inclusive Prosthetic Components labeled STERILE should be stored in a dry location at room temperature, in their original packaging. Users are advised to visually inspect packaging to ensure seals and contents are intact prior to use. Please refer to the individual product label for all relevant product information and cautions. Sterile products are intended for single-use only, prior to the expiration date. Do not use sterile products if the packaging has been compromised or previously opened. Do not resterilize.

Products labeled NON-STERILE should be cleaned and sterilized according to a validated method prior to use in the oral environment.

■ Sterility

Healing abutments and Multi-Unit Abutments are shipped sterile. They should not be re-sterilized. They are for single use only, prior to the expiration date.

Non-sterile abutments and screws must be cleaned and sterilized prior to clinical use, according to a validated method.

- **Cleaning:** Prepare cleaning solution using 5 mL of dish soap per gallon of tap water. Fully immerse the devices in solution and scrub them with a soft-bristle brush. Remove the components and rinse them under running tap water. Dry the devices with a clean, lint-free cloth.

The recommended sterilization process is based on the ANSI/AAMI/ISO 17665-1 and ANSI/AAMI ST79 guidelines, as follows:

- **Sterilization of Inclusive Titanium and Inclusive Esthetic Abutments:** Gravity-fed sterilizers: Autoclave in sterilization pouch for 30 minutes at 121°C (250°F). Devices are to be used immediately after sterilization.
- **Sterilization of Inclusive Hybrid Abutments:** Gravity-fed sterilizers: Autoclave in sterilization pouch for 15 minutes at 132°C (270°F). Allow sterilized components to dry for at least 30 minutes.

NOTE: The validated procedures require the use of FDA-cleared sterilization trays, wraps, biological indicators, chemical indicators, and other sterilization accessories labeled for the sterilization cycle recommended. The healthcare facility should monitor the sterilizer for the facility according to an FDA-recognized sterility assurance standard such as ANSI/AAMI ST79.

■ MRI

The Inclusive abutments have not been evaluated for safety and compatibility in the magnetic resonance (MR) environment. They have not been tested for heating, migration, or image artifact in the MR environment. The safety of Inclusive abutments in the MR environment is unknown. Scanning a patient who has the device may result in patient injury.

■ Dental Implant Compatibility and Recommended Torque Values

Inclusive Prosthetic Components manufactured by PrismaMat Dentalcraft are generally compatible with the implant systems listed in the table below. The availability of a particular type of prosthetic component varies by implant system, and may be limited by geographical territory. The platform-specific compatibility of each component is indicated on the individual product label. For a complete product listing, please refer to the *Inclusive Prosthetic Components Product Catalog*, or contact an Inclusive sales representative.

Implant System and Platform Sizes	Manufacturer's Recommended Torque (Ncm)		
	Titanium Abutment / Screw	Multi-Unit Abutment / Screw	Multi-Unit Prosthetic Screw
Biomet 3i™ Certain® 3.4 mm, 4.1 mm, 5.0 mm, 6.0 mm	20 Ncm	20 Ncm	15 Ncm
Camlog® Screw-Line 3.3 mm, 3.8 mm, 4.3 mm, 5.0 mm, 6.0 mm	20 Ncm	20 Ncm	15 Ncm
Dentsply Implants Ankylos® A (3.5 mm), B (4.5 mm), C (5.5 mm), D (7.0 mm)	25 Ncm	25 Ncm	15 Ncm
Dentsply Implants Astra Tech Implant System® X-Small (2.5 mm) Small (2.9 mm) Large (3.9 mm)	15 Ncm	15 Ncm	15 Ncm
	20 Ncm	20 Ncm	15 Ncm
	25 Ncm	25 Ncm	15 Ncm
Dentsply Implants Astra Tech Implant System® EV 3.0S, 3.6S, 4.2S, 4.2C, 4.8S, 4.8C, 5.4S	25 Ncm	25 Ncm	25 Ncm
Hiossen® HG System Mini (3.5 mm) Standard (4.0 mm)	20 Ncm	20 Ncm	15 Ncm
	30 Ncm	30 Ncm	15 Ncm
MegaGen AnyRidge® Implant System: 3.5 mm, 5.0 mm	30 Ncm	30 Ncm	15 Ncm
Nobel Biocare Brånemark System® RP (3.75 mm)	35 Ncm	35 Ncm	15 Ncm
Nobel Biocare NobelActive® NP (3.5 mm), RP (4.3 mm)	35 Ncm	35 Ncm	15 Ncm
Nobel Biocare NobelReplace®	35 Ncm	35 Ncm	15 Ncm

NP (3.5 mm), RP (4.3 mm), WP (5.0 mm), 6.0 (6.0 mm)			
Straumann® Bone Level NC (3.3 mm), RC (4.1 mm)	35 Ncm	35 Ncm	15 Ncm
Straumann® Tissue Level NN (3.5 mm), RN synOcta® (4.8 mm)*, WN synOcta® (6.5 mm) *excluding 3.3 mmD Tissue Level implants	35 Ncm	35 Ncm	15 Ncm
Zimmer Dental Screw-Vent® 3.5 mm, 4.5 mm, 5.7 mm	30 Ncm	30 Ncm	15 Ncm

NOTE: Any screw-retained prosthetic component not listed in the table above should be hand-tightened only. These values are accurate as of the time of printing. Always verify torque values with the implant manufacturer.

ANKYLOS is a registered trademark of DeguDent GmbH. AnyRidge is a registered trademark of MegaGen Implant Co., LTD. ASTRA TECH Implants is a registered trademark of Dentsply IH AB LLC. BIOMET 3i is a trademark of BIOMET 3i, LLC. Brånemark System is a registered trademark of Nobel Biocare Services AG. CAMLOG is a registered trademark of CAMLOG Biotechnologies GmbH. Certain is a registered trademark of BIOMET 3i, LLC. HIOSSSEN is a registered trademark of OSSTEM IMPLANT Co., LTD. NobelActive is a registered trademark of Nobel Biocare Services AG. NobelReplace is a registered trademark of Nobel Biocare Services AG. Screw-Vent is a registered trademark of Zimmer Dental Inc. Straumann is a registered trademark of Straumann Holding AG. synOcta is a registered trademark of Straumann Holding AG.

■ INSTRUCTIONS FOR USE — INCLUSIVE® ABUTMENTS

Inclusive Abutments are prefabricated, screw-retained intraoral abutments intended to be connected directly to an endosseous implant for retention of a cemented dental prosthesis. They may be indicated for single- and multiple-tooth restorations. Titanium abutments are machined from titanium alloy and attached to the implant fixture with an abutment screw compatible with the restorative instrumentation of the specified implant system.

Capture Implant Placement

Take an implant-level impression utilizing the preferred technique (direct, indirect, or intraoral scan). Submit the impression to the laboratory.

CAD/CAM Preparation

Laboratory — Design the Restoration

- 1) Create a soft tissue model from an implant-level impression.
- 2) Select the appropriate laboratory scanning abutment to capture the implant angulation, position, and abutment connection orientation. Follow manufacturer instructions to obtain all necessary scans to construct an accurate, complete 3-D model.
- 3) Design the abutment according to the patient's clinical needs, taking care to ensure adequate support for the eventual restoration, including appropriate interproximal and occlusal space. Produce a digital design file.
- 4) Send the digital design file to a milling center to manufacture the patient-specific implant abutment.

Milling Center — Fabricate the Restoration

- 1) Select the appropriate Inclusive® Abutment Blank based on the system and platform size of the implant seated in the patient's mouth.
- 2) Fabricate the restoration using CAD/CAM techniques. Veneer as necessary. If a screw-retained hybrid restoration is indicated, fabricate the superstructure (i.e., zirconia coping or crown) and lute it to the titanium abutment. The superstructure is to be bonded to the titanium abutment using MonoCem® Self-Adhesive Resin Cement (Shofu Dental Corporation; San Marcos, Calif.).

Non-CAD/CAM Preparation

Laboratory — Fabricate the Restoration

- 1) Create a soft tissue model from an implant-level impression.
- 2) Select the appropriate Inclusive Abutment based on the system, platform size, location, and occlusal clearance of the implant seated in the patient's mouth.
- 3) Seat the abutment completely into the implant analog on the working model, making sure that the anti-rotational features of the connection interface are fully engaged and the contours of the emergence profile (if applicable) are esthetically oriented.
- 4) Insert the appropriate compatible abutment screw into the abutment's screw access hole and hand-tighten using the appropriate driver.
- 5) Fabricate the restoration using conventional casting techniques. Veneer as necessary. If a screw-retained hybrid restoration is indicated, lute the ceramic crown to the titanium abutment. The ceramic crown is to be bonded to the titanium abutment using MonoCem Self-Adhesive Resin Cement (Shofu Dental Corporation; San Marcos, Calif.).

Manual Adjustment

NOTE: Due to the high thermal conductivity of titanium, titanium abutments should not be modified in the oral cavity. Any necessary modifications should be made extraorally.

- 1) Seat the abutment completely into an implant analog retained by an analog holder or the implant analog captured in the working model, making sure that the anti-rotational features of the connection interface are fully engaged and the contours of the emergence profile (if applicable) are esthetically oriented.
- 2) Insert the appropriate compatible abutment screw into the abutment's screw access hole and hand-tighten using the appropriate driver.
- 3) Using a fine-diamond or carbide bur, modify the abutment as needed.
- 4) With a silicone-based rubber wheel or point, refine the abutment along the margins.

Deliver the Final Restoration

- 1) Seat the titanium abutment or screw-retained hybrid restoration completely into the implant, making sure that the anti-rotational features of the connection interface are fully engaged and the contours of the sculpted emergence profile are esthetically oriented.
- 2) Insert the appropriate compatible abutment screw into the screw access hole and hand-tighten using the appropriate driver. It is strongly recommended that a radiograph of the connection site be taken to confirm complete seating of the abutment or hybrid restoration before proceeding.
- 3) Using the appropriate driver in conjunction with a properly metered torque wrench, tighten the abutment or hybrid restoration to the implant manufacturer's recommended torque value.
- 4) Fill the screw access hole with Teflon tape, gutta-percha, or other suitable material.

5) If the restoration is of a screw-retained hybrid design, cover the screw access hole with flowable composite, and cure. Otherwise, follow applicable cementation procedures to affix the definitive restoration to the abutment.

■ **INSTRUCTIONS FOR USE — INCLUSIVE® TITANIUM ESTHETIC ABUTMENTS**

Inclusive® Titanium Esthetic Abutments are prefabricated, screw-retained intraoral abutments intended to be connected directly to an endosseous implant for retention of a cemented dental prosthesis. They may be indicated for single- and multiple-tooth restorations. Titanium esthetic abutments are precisely machined from titanium alloy and attached to the implant fixture with a titanium screw. Unlike the circular emergence profile of standard stock abutments, esthetic abutments are manufactured with a tapered emergence profile for more natural-looking contouring of the soft tissue at the implant site. Each esthetic abutment is specific to the restorative platform of the implant, and anatomically designed for the connection site's region on the ridge (anterior or posterior). In addition to the standard, straight abutment body, angled abutment bodies, produced with a 15 degree slope of one hemisphere to compensate for an undesirable path of insertion resulting from excessive implant angulation, are available.

Capture Implant Placement

Take an implant-level impression utilizing the preferred technique (direct, indirect, or intraoral scan). Submit the impression to the laboratory.

Laboratory — Fabricate the Restoration

- 1) Follow pouring procedures for the appropriate die stone to produce a working model and articulate with a bite registration.
- 2) Select the appropriate Inclusive Titanium Esthetic Abutment based on the system, platform size, location, angulation, and occlusal clearance of the implant seated in the patient's mouth.
- 3) Seat the abutment completely into the implant analog on the working model, making sure that the anti-rotational features of the connection interface are fully engaged and the contours of the sculpted emergence profile are esthetically oriented. For angled abutments, the tapered side should be oriented nearest vertical along the same plane as the implant.
- 4) Insert an Inclusive Titanium Screw into the abutment's screw access hole and hand-tighten using the appropriate driver.
- 5) Fabricate the restoration using conventional casting or CAD/CAM techniques. Veneer as necessary. If a screw-retained hybrid restoration is indicated, lute the ceramic crown to the titanium abutment.

Deliver the Final Restoration

- 1) Seat the titanium esthetic abutment or screw-retained hybrid restoration completely into the implant, making sure that the anti-rotational features of the connection interface are fully engaged and the contours of the sculpted emergence profile are esthetically oriented. For angled abutments, the tapered side should be oriented nearest vertical along the same plane as the implant.
- 2) Insert an Inclusive Titanium Screw into the screw access hole and hand-tighten using the appropriate driver. It is strongly recommended that a radiograph of the connection site be taken to confirm complete seating of the abutment or hybrid restoration before proceeding.
- 3) Using the appropriate driver in conjunction with a properly metered torque wrench, tighten the abutment or hybrid restoration to the implant manufacturer's recommended torque value.
- 4) Fill the screw access hole with Teflon tape, gutta-percha, or other suitable material.
- 5) Follow applicable cementation procedures to affix the definitive restoration to the abutment. Or, if the restoration is of a screw-retained hybrid design, cover the screw access hole with flowable composite, and cure.

■ **INSTRUCTIONS FOR USE — INCLUSIVE® MULTI-UNIT ABUTMENTS**

Inclusive Multi-Unit Abutments are prefabricated, screw-retained intraoral abutments intended to be connected directly to endosseous implants in partially or fully edentulous patients for the retention of cast or milled bar overdentures. Multi-unit abutments are machined from titanium alloy and are available with a variety of collar heights to achieve optimal emergence from shallow or deep gingival wells. Each Inclusive Multi-Unit Abutment is delivered sterile, with a carrier color-coded to indicate the restorative platform of the implant.

Straight Multi-Unit Abutments lack any anti-rotational features at the implant-abutment interface. The apical portion of a straight Multi-Unit Abutment is threaded for integration with the internal cavity of the implant. For abutment delivery, the occlusal surface features a male hex head compatible with the Inclusive Multi-Unit Driver.

Angled Multi-Unit Abutments of 17 degrees or 30 degrees enable clinicians to compensate for the divergence of implants or to otherwise accommodate an angled path of insertion. Angled Multi-Unit Abutments feature an anti-rotational connection interface specific to the matching implant platform, and are attached to the implant fixture with an angled Multi-Unit Abutment screw compatible with the restorative instrumentation of the specified implant system. Both straight and angled Multi-Unit Abutments feature a female connection port at the coronal apex, to allow for the attachment of a screw-retained or fixed-removable dental prosthesis with a multi-unit restorative screw (Inclusive Prosthetic Screw).

The axial tilt of an Inclusive Angled Multi-Unit Abutment (angular divergence from path of insertion) is designed and manufactured to lie along a *plane* of the implant connection geometry, as opposed to a corner or junction. To maximize the angle-correcting attributes of the Multi-Unit Abutment, be sure to rotate the implant upon final seating so that one side of the internal connection geometry (flat or lobe) is oriented to serve as the base of angulation, in accordance with the restorative treatment plan.

Place the Multi-Unit Abutment

- 1) Select the appropriate Inclusive Multi-Unit Abutment based on platform size, implant angle, and depth of the soft-tissue well.
- 2) Retrieve the abutment from its packaging. To maintain the sterility of the Multi-Unit Abutment, be careful to handle only by the carrier.
- 3) (a) *For Straight Abutments:* Using the carrier, seat the abutment into the implant and hand-tighten. Remove the carrier by pulling the apex of the carrier toward the facial. (b) *For Angled Abutments:* Using the carrier, seat the abutment into the implant until the anti-rotational features of the connection interface are engaged. Lift and rotate as necessary to orient the angle in the required direction. Hand-tighten the Inclusive Angled Multi-Unit Abutment Screw using the appropriate driver. Twist the carrier counterclockwise to unscrew the carrier from the abutment.

NOTE: It is strongly recommended that a radiograph of the connection site be taken to confirm complete seating of the abutment before proceeding.

- 4) Using the appropriate driver in conjunction with a properly metered torque wrench, tighten the Multi-Unit Abutment or angled Multi-Unit Abutment screw to the implant manufacturer's recommended torque value.

Passive Temporization of Multi-Unit Abutments

- 1) If delayed loading is desired, cover each Inclusive Multi-Unit Abutment with an Inclusive Multi-Unit Temporary Healing Cap and hand-tighten with the Inclusive Prosthetic Screw provided, using the appropriate driver. Do not overtighten.

- 2) Using the patient's existing denture or other prosthesis, relieve the area directly above the placement of each temporary healing cap until the denture rests on the ridge.
- 3) Follow procedures to reline the denture over the temporary healing caps, using tissue conditioner or soft relining material. The temporized denture can be used during healing.

NOTE: For a temporization technique involving loading, please refer to the *Inclusive Prosthetic Components Restorative Manual*.

Capture Multi-Unit Abutment Placement

Take an abutment-level impression utilizing the preferred technique (direct, indirect, or intraoral scan). Submit the impression to the laboratory for the fabrication of a working cast and verification index.

Denture Protocol

Follow appropriate denture protocol in accordance with the patient-specific treatment plan. When trying in the various setups (e.g., verification index, occlusal rim, wax setup, retention bar), hand-tighten to the Multi-Unit Abutments with prosthetic screws, using the appropriate driver. Start from the distal and move forward, alternating between sides of the ridge. Always confirm complete, passive seating, modifying the setup as needed.

Deliver the Final Restoration

- 1) Remove any temporary prosthesis.
- 2) Confirm that each Multi-Unit Abutment is tightened to the implant manufacturer's recommended torque value.
- 3) Align the prosthesis onto the abutments. Beginning with the midmost screw access channel, hand-tighten an Inclusive Prosthetic Screw into the Multi-Unit Abutment. Repeat for each abutment, working outward and alternating left to right.
- 4) Confirm appropriate seating. With the same middle-out, left-to-right technique, tighten each prosthetic screw to the recommended torque value.
- 5) Check comfort and occlusion, and make any necessary adjustments.
- 6) Place Teflon tape, gutta-percha, or other suitable material.
- 7) Cover each screw access hole with flowable composite and cure.

SYMBOLS GLOSSARY

Symbol	Symbol Ref. No.	Symbol Title	Designation No.	Explanatory Text
	5.2.4	Sterile with Gamma Radiation	EN ISO 15223-1	This symbol indicates that this device has been sterilized using irradiation.
	5.2.8	Do Not Use if Package is Damaged	EN ISO 15223-1	This symbol indicates that this device should not be used if the package has been damaged or opened.
	5.2.7	Non-Sterile	EN ISO 15223-1	This device has not been subjected to a sterilization process.
	5.4.2	Do not Re-use	EN ISO 15223-1	This device is intended for one use, or for use on a single patient during a single procedure.
	5.2.6	Do not Resterilize	EN ISO 15223-1	This symbol indicates that this device is not to be resterilized.
	5.1.4	Use-by Date	EN ISO 15223-1	This symbol indicates the date (YYYY-MM-DD) after which this device is not to be used.
	Sec. 801.109(b)(1)	By Prescription Only	21 CFR Part 801	Caution: Federal law restricts this device to sale by, or on the order of, a licensed dentist or physician.
	5.1.6	Catalog Number	EN ISO 15223-1	This symbol indicates PrismaTik Dentalcraft's catalog number so that this device can be identified.
	5.1.5	Lot/Batch Number	EN ISO 15223-1	This symbol indicates PrismaTik Dentalcraft's lot/batch number so that the lot/batch of this device can be identified.
	5.4.3	Consult Instructions For Use	EN ISO 15223-1	This symbol indicates the need of the user to consult the instructions for use.
	5.1.1	Manufacturer Date of Manufacture (YYYY-MM-DD)	EN ISO 15223-1	This symbol indicates the manufacturer and the date of manufacture of this device.

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Prismatic Dentalcraft, Inc.